

eBIC for Industry

Flexible cryo-EM services to suit all levels


The electron Bio-Imaging Centre (eBIC) at Diamond Light Source offers bespoke services to industry clients wanting to learn or de-risk the adoption of cryo-EM in their research projects. By offering flexible access to world-class cryo-EM facilities, including dedicated Glacios and Krios microscopes, and highly-experienced scientists, eBIC for Industry lowers the barrier to entry for a technique traditionally requiring large capital and training investments.



Situated within the integrated synchrotron environment, our services represent a unique offer to industry, with a cross-over between structural techniques like MX and SAXS as well as dedicated experts for industry, and streamlined support services.


On-demand access to world-class facilities and expertise

- Industry-dedicated instruments and laboratories
- Full scientific consultancy and 1:1 support from experienced cryo-EM scientists
- Streamlined user, data, logistics, and legal functions
- Objective assessment of project suitability and progress




State-of-the-art capabilities for data collection

- Screening and preliminary data collection on our Glacios for rapid feedback
- High-throughput and best resolution data collection on our industry-dedicated Krios
- EPU multi-grid data collection and on-the-fly processing



	SCREENING	COLLECTION	
Microscope	Glacios™	Krios™ G3i	
Detector	TFS Falcon4i	TFS Falcon4i + SelectrisX	
Mode	Counting	Counting	
Acquisition software	EPU (3.3)	EPU (3.4)	
Image Size (pixels)	4096 x 4096	4096 x 4096	
Magnification (kx)	120	105	130
Pixel size (Å)	0.85	1.177	0.925
Dose rate (e-/px/s)	8.8	15.0	6.0
Total Dose (e-/Å²)	49	38	37
Exposure time (s)	4	3.5	5.23
EER Internal Frames	1224	1080	1620
Acquisition Rate (per hour) R1.2 type	400	566	471
Movie Size (MB)	569	650	520
Movie Format	EER	EER	
Average Dataset Size (TB)	6	4.5 - 6.5	

Data Quality vs. Quantity for eB4i Collection Configurations



Flexible and scalable access modes


- 1 Rental:** come and use our labs and microscopes in person
- 2 Remote:** send grids and log into the microscopes to collect the data yourself
- 3 Mail-in:** send your grids and our experts will collect data for you
- 4 Full-service:** we freeze, screen, collect, and solve your structures



Services tailored to meet your requirements

Full support for companies adopting cryo-EM

- From start-ups without labs to existing microscope owners
- Lower upfront costs and de-risking of EM adoption, with scalable access
- Clients can participate in all aspects of the experiment
- Industrial scientists can enter or leave the pipeline at any stage




Applied training at any level

- Bespoke, hands-on, practical training in sample preparation and data collection
- Progress real samples whilst learning EM
- Attend in-person or remotely



Rapid structural insights for challenging proteins

- SPA structures of ligand/fragment-bound proteins
- Macromolecular complexes and flexibility analysis
- IHRSR of amyloid structures
- Rational vaccine design
- Antibody epitope mapping
- Drug delivery system and gene-therapy vector characterisation



For further information
Diamond Industrial Liaison Team

+44 1235 778797
industry@diamond.ac.uk

diamond.ac.uk/industry
@DiamondILO