

Re-Window an XRF map

What this does:

To produce the XRF maps during an experiment you specify the elements you're interested in and the software sums N data channels around the XRF peak in your data.

If you've forgotten to include a peak or want to adjust the number of channels used in generating a map you can re-process your data.

Note that windowing does not deconvolute overlapping XRF peaks so should be treated with caution. For example in the case of Fe and Co the Fe K β transition occurs close to the energy of the Co K α , windowing for Co.


Open Dawn

In Linux - Open a terminal then type:

```
module load dawn
dawn&
```


Alternatively, you can open Dawn double-clicking the corresponding "script-shorcut" present in the desktop if you have previously set-up the desktop configuration

(typing in Terminal: "`cd /dls_sw/i14/scripts/user_setup/configure_desktop/`" — "Enter" — Then typing: "`./setup_user.sh`")




In Windows - launch dawn from your start menu


Open the Processing perspective:



OR




Select the processing perspective:



Load your data set:


In the top-left side of the program, select the "File Viewer" tab. Scroll down and select your experimental folder (e.g. cm22977-4) in the corresponding path (e.g. /dls/14/data/2019).

Then, go to the "scan" folder and click the file which need reprocessing. The "File Filter" option can be used to find a specific scan easier, by typing the scan number and extension (*59397.nxs).



Select the file you want and double-click it or drag it into the Data Slice View window. A dialog will open with a drop-down menu where you have to select the "/entry/xsp3_detector/data" option.


This menu should have selected the detector data and X and Y values correctly.



However - **The type of the DIM1 and DIM2 values may need to be "RANGE"** in some cases, as shown in the screenshot below.


Set up data for processing

Select dataset, axes, whether to process as images [2D] or lines [1D] and which dimensions of the array are the data dimensions




Press Finish

In the processing window click on the start and type XRF - this should select the 'XRF Elemental Maps from ROIs' tool and press return.



Press the live setup option at the top-right corner:



Select the elements and lines you want and set the ROI width as needed. Beforehand, remember to add the value '20' as ROI width (# channels) at the bottom right, and to select the excitation energy used so only the excited elemental transitions can be selected.

Region of interest based elemental maps of 2D XRF datasets

Periodic Table Analyze existing spectrum

Selected XRF linegroups

Fe-Kα

Remove linegroup

Min excitation energy (keV) 1.0

Max excitation energy (keV) 50.0

ROI width (# channels) 10

Cancel Finish

You can examine the ROI width and raw MCA in the Analyze Tab

Region of interest based elemental maps of 2D XRF datasets

Periodic Table Analyze existing spectrum

Fe - 26

Intensity

Energy (keV)

— Fe — Sum peaks

No escape peaks Show sum peaks Energy unit is 1 keV

Selected XRF linegroups

Fe-Kα

Remove linegroup

Min excitation energy (keV) 1.0

Max excitation energy (keV) 50.0

ROI width (# channels) 10

Use current dataset
 Use sum dataset
 Use external dataset


Kα
 Kβ
 Lα
 Lβ
 Mα

Cancel Finish

Press finish once you have set it up.

Then press the green play button - "Process all files" - on the "Data Slice View" tab. This will subsequently open a dialog before processing the selected file.

Ensure that the copy directory is the "processing" folder of your visit (NOT in the 'processed' as shown in the image below, since this saving files to this folder is restricted). Select the type of processing data required, 'processed data only' would generate a lighter file with the new data only, but you can either 'link the original data' in the new file or create a new file with all the information (last option available, which will take longer because of the heavier file to be created).



Subsequently click OK.

Processing more files or examining your processed data:

If you want to run lots of files you can drag them into the "Data Slice View" tab.

Alternatively - once you have setup a processing chain and run it one - double clicking on any file in processing perspective will automatically process that file in the same manner.

The processing perspective is to process files it is not for viewing the results. You would need to look at your data in the "Mapping" perspective.

File Tools Window Help

Quick Access DataVis Mapping Processing

Project File View

/2019/cm22977-4/scan

All Folders

- i12
- i13
- i13-1
- i14
 - snapshots
 - data
 - 2015
 - 2016
 - 2017
 - 2018
 - 2019
 - cm22977-1
 - cm22977-2
 - cm22977-3
 - cm22977-4
 - .ispyb
 - .workspace
 - processed
 - processing
 - scan
 - spool

Contents of /dls/i14/data/2019/cm22977-4/scan

File Filter *59397.nxs

Name	Size	Type
i14-59397.nxs	279.6 KiB	N

1 object(s) selected 286.35

Data Slice View


Process all files

- i14-59397.nxs
- i14-59398.nxs
- i14-59399.nxs
- i14-59400.nxs
- i14-59401.nxs
- i14-59402.nxs
- i14-59403.nxs
- i14-59404.nxs
- i14-59405.nxs
- i14-59406.nxs
- i14-59407.nxs
- i14-59408.nxs
- i14-59409.nxs

Current slice of data: [0,0,4,4096]

Input

detector[0,0,::]



Output

XRF Elemental Maps from ROIs [XRF Analysis]

Name	Value
------	-------

