

XChem lab checklist

- Have your Diamond FedID and password in hand
- Check if there is enough pucks/pins for your experiment and recycle if needed
- Image your plate at the Research Complex
- Set xchempath in your linux bashrc if you are new on the platform
- Run: preparevisit to create the XChem visit directories
- Open SoakDB: /dls/labxchem/data/PROPOSAL/VISIT/processing/lab36/
- Input your protein acronym (this must match what is in your ERA)
- Rank your crystals using TexRank (remember to set to Formulatrix and 1 drop) and export the target list CSV file
- Watch out for X/Y coordinates > +/-1400 in the exported CSV file
- If using plate ≠ SwissCi-3drop, change the plate type in SoakDB
- Export Echo soaking CSV files from SoakDB
- Open compound plate under argon gas
- Transfer compounds for soaking to your crystal drops using the Echo
- Flush compound plate with argon gas before sealing it
- Update the transfer status to “Done” in SoakDB
- Export the harvesting CSV file for the shifter
- Scan the required pucks out at the scan in/out workstation
- Scan your pucks on the Shifter software in the same order you are going to use them and select the first row of the puck list
- Mount your crystals
- Scan your open pucks/pins on the puck scanner station
- Seal your pucks and store in the dewars and write down the cane tape number
- If using plate ≠ SwissCi-3drop, before you import the harvesting CSV file back, be sure you still have the right plate type selected in SoakDB
- Import your harvesting CSV shifter file back into SoakDB
- Click on “update status” in the SoakDB ISpyB section
- Add beamline visit, unit cell dimensions and data collection parameters for UDC in SoakDB as necessary
- Click on “Export to ISpyB” in the SoakDB ISpyB section
- Open Putty and connect to `ssh.diamond.ac.uk` (click on “yes” if a pop-up window asking to add the RSA key appears)
- Login with your FedID and password and go to your visit ispyb directory: /dls/labxchem/data/PROPOSAL/VISIT/processing/lab36/ispyb/
- Run the script: `csv2ispyb your-exported-ispyb-csv-file.csv`
- Log onto ISpyB and check that a new shipment with your pucks has been created
- Go to SoakDB “Pucks” tab, click on “Update/Save changes” then on “Yes”. Update the puck table with: the puck location, the data collection parameters, the beamline visit number, your name, your target name and your local contact
- Copy the new/relevant lines into the XChem-Queue Google spreadsheet
- Close your SoakDB
- Scan your unused pucks back in the scan in/out workstation