



I05 – Angle Resolved Photoemission Spectroscopy (ARPES)

Angle-Resolved Photoemission Spectroscopy (ARPES) maps the dispersion of electronic bands near the Fermi level and, in particular, the Fermi surface itself by exciting the bound electrons in a sample with a given photon energy. The momentum parallel to the surface is fully conserved, thus making the method suitable for layered low-dimensional materials. The three-dimensional momentum distribution is also reflected in the photoelectron features thus making the spectroscopy applicable to metallic single crystals, provided that a well-defined clean surface can be prepared in ultra-high vacuum. The minimum samples size is $500 \times 500 \mu\text{m}^2$ given by the light spot ($50 \times 50 \mu\text{m}^2$) and the sphere of confusion of the sample goniometer.

Sample Preparation

Samples have to be prepared *in situ* in ultrahigh vacuum to achieve atomically controlled, clean surfaces. The vacuum system of the high resolution ARPES branch provides the following methods:

- cleavage of layered materials by knocking off a post glued to the top surface
- Ar-ion etching and annealing for polished metal surfaces
- evaporation sources for the preparation of ultrathin films
- direct heating flash for passivated semiconductor surfaces
- port for attaching user-supplied preparation equipment



Beamline Specification

Beamline I05	High Resolution Branch
Photon Energy Range	18 – 240 eV
Energy Resolution	10 meV
Angular Resolution	0.1°
Available Polarisation Modes	Linear Horizontal/Vertical, Circular Left/Right
Smallest Sample Size	$500 \times 500 \mu\text{m}^2$
Controlled Sample Temperature Range	10 – 400 K
Vacuum Conditions	$< 5 \times 10^{-11}$ mbar

For further information please contact the Diamond Industrial Liaison Office on



+44 (0)1235 778797



industry@diamond.ac.uk




www.diamond.ac.uk/industry



@DiamondILO


Carbon-Based Materials

- Graphene
- Molecular electronics
- Carbon nanotubes




New Materials

- ARPES as soon as single crystals are grown




Transition Metals

- Quantitative analysis of electron interactions
- Fermi surfaces
- Renormalisation
- Energy gaps



Surfaces and Interfaces

- Molecular adsorbates
- Ultrathin films
- Stepped surfaces
- Epitaxially grown nano-wires
- Topological insulators



For further information please contact the Diamond Industrial Liaison Office on



+44 (0)1235 778797



industry@diamond.ac.uk



www.diamond.ac.uk/industry



@DiamondILO