LCLS X-ray mirror measurements using a large aperture visible light interferometer

T. J. McCarville, R. Soufli, M. J. Pivovaroff Lawrence Livermore National Laboratory

ACTON 2011 conference



This work performed under the auspices of the U.S. Department of Energy by Lawrence Livermore National Laboratory under Contract DE-AC52-07NA27344 in support of LCLS Project at Stanford Linear Accelerator Center.

This talk describes how we arrived at mirror figure tolerances for LCLS hard X-ray mirrors, and how they where measured



Nominal high energy beam line layout

Mirror type	Tangential Sphere (peak/valley)	Tangential A-sphere (rms)
Flats	< +/- 10 nm	< 2 nm
K/B focus mirror	< +/- 3 nm	< 1 nm

The affect of figure errors on phase/intensity depends on the figure error length scale

• A figure error of length L(0) initiates a phase disturbance of length $L(0)\sin\Theta$

• The disturbance grows in size as $L(z) = L(0)(\sin \Theta) [1 + (z/Ra)^2]^{1/2}$, where the Rayleigh length $Ra = \pi (L(0)\sin \Theta)^2 / (4\lambda)$



- (1) rapidly expand and intermingle
- (2) gradually evolve in both phase and intensity
- (3) propagate as a phase-only disturbance
- Most figure errors of interest lead to gradual evolution
 - examined by solving the Fresnel equation



L(0)

Θ

L(z)

A tolerance for spherical error in flats and focus mirrors can be arrived at analytically – using peak intensity at focus as the metric

• An Optical Path Difference (OPD) results when mirror spherical curvature differs in sagittal and tangential planes by an amount δ , peak/valley



• For astigmatism, the relative peak intensity I at focus scales as (Born & Wolf, sec. 9.3): $I > 1 - \frac{\pi^2}{6\lambda^2} \sin^2(2\Theta\delta)$

Maintaining < 80% peak intensity limits the spherical error in flats to < 6 nm (peak/valley)
< 2 nm relative error between focus mirrors



Once the spherical error is known, it can be reduced to a focus shift Δf by translating K/B focusing optics relative to one another



 $\Delta f/f = \theta/NA$ • In terms of flat and focus mirror geometry, $\Delta f = 16f^2(\delta/L)_{flat} / (L\Theta)_{focus}$ = 7 mm shift per nm of spherical error in an LCLS flat $0 = 8(\delta/L)_{flat}$ • the 10 nm of sphere allowed in the flats can be corrected by a shifting a K/B mirror 70 mm

Measuring sphere accurately enables a focus shift to be pre-planned prior to mirror installation

Tolerances for a-sphere in LCLS flats were specified using the Maréchal criterion

• The criterion describes the reduction in peak intensity at focus due to a random figure error distribution of rms amplitude δ_{rms}

$$I < 1 - \left(\frac{2\pi}{\lambda}\right)^2 (2\delta_{rms}\sin\Theta)^2$$

• The error must be < 3 nm rms to achieve > 80% of the diffraction limited intensity at focus

• A smaller number applies to focus mirrors, due to larger graze angle



• Measuring these tolerances with a visible light interferometer requires:

- Measurement noise << 1 nm (must be a phasing interferometer)
- Normal incidence, using only the essential optics (no turning flats, etc.)
- Careful attention to mirror distortion while mounting
- Absolute calibration to < 1 nm with a three-flat test

A three flat test was used to calibrate a 300 mm diameter interferometer transmission flat to < 0.5 nm (rms) along the horizontal axis

• The figure of transmission flats T1, T2, and reflection flat R were solved along a horizontal line using the geometries shown



• The solution error was sampled using two additional independent configurations





The a-sphere figure error varies from 1 to 2.6 nm (rms) for four 450 mm long flats

- The flats are fabricated by Zeiss with 150 nm of concave spherical curvature, and < 2 nm rms a-sphere
 - they are then coated by Regina Soufli at LLNL, reducing the curvature by about 10 nm
 - the sphere is mechanically reduced to < 10 nm peak/valley during mounting, in front of an interferometer
 - a-sphere is monitored throughout to make sure it is not increased my mounting forces





- Zeiss measures the a-sphere independently prior shipping
 - Zeiss & LLNL measurements are in good agreement



Focus mirrors with high fringe density can be accurately measured by calibrating the "trace back" error

- Pixel density limits the maximum curvature that can measured
 - Nyquist limit requires > 3 pixels/fringe
 Limits measurements to < 25 microns (peak/valley)
 across the aperture
- At high fringe density, sample rays and reference rays travel different paths through the interferometer
 - a measurement error is observed when a test flat is tilted off null fringe





• The correction that would be applied to the f=8.2 m K/B mirror is shown

- note a slight miss-alignment in the imaging optics
- the correction is small, but significant compared to
- 1 nm tolerance for a-sphere



Conclusions and Future Efforts

• A visible light interferometer was used to accurately measure 450 mm long flats for LCLS

- bend mirrors flat to < 10 nm peak/valley

- measured a-spheric figure ranging from 1 - 2.6 nm rms for four flats

• Focus mirrors can be measured accurately by applying a correction for systematic errors of the interferometer imaging optics

- the correction is significant compared to focus mirror figure tolerances

- The interferometers 300 mm aperture requires file stitching, which limits measurement accuracy
 - plans are underway to calibrate a 600 mm instrument
 - 900 mm aperture instruments are becoming available

Acknowledgements

The authors wish to thank Peter Stefan and Jacek Krzywinski at Stanford Linear Accelerator Center for their many contributions and encouragement.