

## CCP4 I2 Getting Started


To open the CCP4 I2 graphical interface, open a terminal and type the following command

```
$ ccp4i2
```

This will open the CCP4 I2 welcome screen.

To create a project, either click the link labelled “Start a new crystallographic project” or select “Projects -> New project” from the application menu. This will open a dialog window asking for the project name. Enter a name for your project. (N.B. pressing “Enter” will select “Cancel” so you must use a mouse here).

By default, the project will be created in your home directory. To save the project in your visit directory, choose “Select directory” and then navigate your visit directory. This is necessary because your home directory has a limited memory quota.



Once the project has been created, the main CCP4 I2 window will then open.

